

Essential Angle Feedback



- Designed to improve joint control in various positions
- Assists with functional training by providing audio and visual feedback on joint movement e.g. knee hyperextension
- Portable, easy to set up, simple and efficient to use daily during rehabilitation
- Angle sensor can be placed over the knee, hip and spine and the elbow using straps
- Angle sensor is connected via a cable to a simple control unit
- Therapist can establish an upper and a lower limit
- Client is given real time audio feedback via a speaker or headphone when the joint angle is outside these limits
- Real time visual feedback is provided by a 10 coloured LED display
- Designed and made in Australia

Optional Extra: DataQ hardware and software to be used with Essential Angle Sensor Unit

DataQ hardware and software have been pre-configured to accept, record and save data from either the Essential Angle or Essential Seating control units. The control unit is connected via a cable to the DataQ box. Once the DataQ box is connected to a computer, the data can be displayed as a graph, recorded and saved to Excel.

